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Reader Display	Interpretation
Flu A: Positive Flu B: Negative Procedural Control: Valid	Positive Test for Flu A (influenza A antigen present)
Flu A: Negative Flu B: Positive Procedural Control: Valid	Positive Test for Flu B (influenza B antigen present)
Flu A: Positive Flu B: Positive Procedural Control: Valid	Positive Test for both Flu A and Flu B (influenza A and B antigen present)
Flu A: Negative Flu B: Negative Procedural Control: Valid	Negative Test for Flu A and Flu B (no antigen detected)
Flu A: Invalid Flu B: Invalid Procedural Control: Invalid	Result Invalid

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- 2) Lee CK, Cho CH, Woo MK, Nyeck AE, Lim CS, Kim WJ. Evaluation of Sofia fluorescent immunoassay analyzer for influenza A/B virus. *Journal of Clinical Virology* 2012; 55(3): 239-243.
- 3) Rath B, Teif F, Obermeier P, Tuerk E, Karsch K, Muehlhans S, Adamou E, et al. Early detection of influenza A and B infection in infants and children using conventional and fluorescence-based rapid testing. *Journal of Clinical Virology* 2012; 55(4): 329-323.

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I. Device Description:

The Sofia Influenza A+B FIA employs immunofluorescence technology that is used with the Sofia Analyzer to detect influenza virus nucleoproteins.

The Sofia Influenza A+B FIA is a lateral-flow immunoassay that uses monoclonal antibodies that are specific for influenza antigens and have no known cross-reactivity to normal flora or other known respiratory pathogens.

Nasal swab, nasopharyngeal swab, and nasopharyngeal aspirate/wash specimens are used for this test. The patient specimen is placed in the Reagent Tube, during which time the virus particles in the specimen are disrupted, exposing internal viral nucleoproteins. After disruption, the specimen is dispensed into the cassette sample well. From the sample well, the specimen migrates through a test strip containing various unique chemical environments. If influenza viral antigen is present, they will be trapped in a specific location.

The Sofia Analyzer will scan the test strip and measure the fluorescent signal by processing the results using method-specific algorithms. The Sofia Analyzer will display the test results (Positive, Negative, or Invalid) on the screen.

L. Test Principle:

The Sofia Influenza A+B FIA is an immunoassay that uses a sandwich design to detect and differentiate Influenza A and Influenza B in patient specimens. This test uses a lateral flow design with location-dependent lines and zones. The Sofia Analyzer scans the test strip and displays the results (Positive, Negative or Invalid) after utilizing the method-specific algorithms. Basically Sofia Influenza A+B FIA involves the disruption of Influenza A and B viral antigens and exposure of internal nucleoproteins by placing the patient specimen in the reagent tube. After disruption the specimen is dispensed into the Cassette sample well after which the sample is drawn by capillary action into and through the label pad, through the nitrocellulose strip and into the adsorbent pad. Within the label pad, the specimen comes into contact with antibodies that have been coupled to the europium chelate-impregnated microbeads. During this interaction, the independent subpopulations of beads containing bound anti-influenza A or anti-influenza B monoclonal antibodies bind corresponding influenza A or B nucleoprotein antigens that are present in the specimen. The bead-coupled antigen-antibody complexes then begin to flow through the test strip. Some of these beads bind non-specifically to the NC line. Most migrate on and, if influenza A or B antigen are present in the sample, they will be subsequently captured on the surface of the nitrocellulose by the respective location-fixed, analyte-specific anti-influenza A or anti-influenza B capture antibodies. The flow and capture of the fluorescent microbeads coated with influenza A and/or B nucleoprotein antigens allows the accumulation of a fluorescent signal at specific analyte line locations on the test strip. The conjugates that do not bind to the Negative Control line and Test Lines continue to flow with the remaining specimen and soon encounter the Reference line that is comprised of goat anti-mouse Ig. The fluorescent signal generated at this line serves as the location marker to direct the Sofia Analyzer to the other precise locations on the Nitrocellulose that are to be scanned. The remaining sample then flows into the Procedural Control Zone that is also scanned to confirm that adequate flow of the sample has occurred. The Sample with any remaining conjugate then flows on into the Adsorbent Pad. This process takes approximately fifteen minutes, at which time the Sofia Analyzer scans the strip, measuring the fluorescent signals across the strip's length and performs calculations and reports the test results. The assay's sensitivity is derived from the use of a unique polystyrene microbead that has been dyed with a chelate of europium that is temperature stable, resistant to bleaching in room light and yields a very efficient conversion of the UV energy from 365 nm to a wavelength of 618 nm.

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1) 미국 질병통제예방센터(Centers for Diseases Control and Prevention, CDC) - Overview of Influenza Testing Methods (last reviewed: 2019.4.4.)

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Method ¹	Type Detected	Acceptable Specimens ²	Test time	CLIA Waived ³
Rapid Influenza Diagnostic Tests ⁴ (antigen detection)	A and B	NP ⁵ swab, aspirate or wash, nasal swab, aspirate or wash, throat swab	<15 min.	Yes/No
Rapid Molecular Assay [influenza viral RNA or nucleic acid detection]	A and B	NP ⁵ swab, nasal swab	15-30 minutes ⁶	Yes/No ⁶
Immunofluorescence, Direct (DFA) or Indirect (IFA) Florescent Antibody Staining [antigen detection]	A and B	NP ⁵ swab or wash, bronchial wash, nasal or endotracheal aspirate	1-4 hours	No
RT-PCR ⁷ (singleplex and multiplex; real-time and other RNA-based) and other molecular assays [influenza viral RNA or nucleic acid detection]	A and B	NP ⁵ swab, throat swab, NP5 or bronchial wash, nasal or endotracheal aspirate, sputum	Varies (1 to 8 hours, varies by the assay)	No
Rapid cell culture (shell vials; cell mixtures; yields live virus)	A and B	NP ⁵ swab, throat swab, NP5 or bronchial wash, nasal or endotracheal aspirate, sputum; (specimens placed in VTM ⁸)	1-3 days	No
Viral tissue cell culture (conventional; yields live virus)	A and B	NP ⁵ swab, throat swab, NP5 or bronchial wash, nasal or endotracheal aspirate, sputum (specimens placed in VTM ⁸)	3-10 days	No

¹ Serologic (antibody detection) testing is not recommended for routine patient diagnosis and cannot inform clinical management. A single acute serum specimen for seasonal influenza serology is uninterpretable and should not be collected. Serological testing for detection of antibodies to seasonal influenza viruses is useful for research studies and requires collection of appropriately timed acute and convalescent serum specimens and testing of paired sera at specialized research or public health laboratories.

² Approved clinical specimens vary by influenza test. Consult the manufacturer's package insert for the approved clinical specimens for each test. Ref: Leland, et al. 2007, Clin Micro Rev 20: 49-78. Approved respiratory specimens vary among FDA cleared influenza assays.

³ Clinical Laboratory Improvement Amendments (CLIA) of 1988.

⁴ **Chromatographic- and/or fluorescence-based lateral flow and membrane-based immunoassays.** Some approved rapid influenza diagnostic assays utilize an analyzer reader device.

⁵ NP = nasopharyngeal

⁶ Rapid molecular assays can provide results in approximately 15-30 minutes.

⁷ Reverse transcription polymerase chain reaction, including FDA-approved test systems, reference laboratory testing using ASR or lab-developed reagents. Some approved molecular assays can produce results in approximately 60-80 minutes.

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Procedure (Manufacturer/Distributor)	Influenza Virus Types Detected	Approved Specimens ²	CLIA Waived ³
3M™ Rapid Detection Flu A+B Test ⁴ (3M)	A and B	NP ⁵ swab/aspirate Nasal wash/aspirate	No
Alere™ Influenza A & B ⁴ (Alere)	A and B	Nasal swab	Yes
BD Veritor™ System for Rapid Detection of Flu A+B CLIA-waived ^{4,8} (Becton Dickinson)	A and B	NP ⁵ swab/nasal swab	Yes
BD Veritor™ System for Rapid Detection of Flu A+B Moderately Complex ^{4,8} (Becton Dickinson)	A and B	NP ⁵ wash/aspirate/swab	No
BinaxNOW® Influenza A&B ⁴ (Alere)	A and B	NP ⁵ swab, Nasal wash/aspirate/swab	Yes
BioSign® Flu A+B ⁴ (Princeton BioMedtech)	A and B	NP ⁵ swab/aspirate/wash, nasal swab	No
Directigen™ EZ Flu A+B ⁴ (Becton-Dickinson)	A and B	NP ⁵ wash/aspirate/swab Throat swab	No
OSOM® Influenza A&B ⁴ (Sekisui Diagnostics)	A and B	Nasal swab	No
QuickVue® Influenza Test ⁶ (Quidel)	A or B	Nasal wash/aspirate/swab	Yes
QuickVue® Influenza A+B Test ⁴ (Quidel)	A and B	NP ⁵ swab Nasal wash/aspirate/swab	Yes
SAS™ FluAlert A&B (SA Scientific)	A and B	Nasal wash/aspirate	No
SAS™ FluAlert A ⁴ (SA Scientific)	A only	Nasal wash/aspirate	Yes
SAS™ FluAlert B ⁵ (SA Scientific)	B only	Nasal wash/aspirate	Yes
Sofia® Influenza A+B ^{4,7,8} (Quidel)	A and B	NP ⁵ aspirate/swab/wash Nasal wash	Yes
TRU FLU® ⁴ (Meridian Bioscience)	A and B	NP ⁵ aspirate/swab Nasal wash	No
XPECT™ Flu A&B ⁴ (Remel/Thermofisher)	A and B	Nasal wash/swab Throat swab	No

1. List may not include all test kits approved by the U.S. Food and Drug Administration. Discontinued tests not included.

2. Approved respiratory specimens according to manufacturer's package insert. Note that test performance may vary if other respiratory specimens are used.

3. Ref: <http://www.cms.gov/Regulations-and-Guidance/Legislation/CLIA/index.html>

4. Distinguishes between influenza A and B virus infections.

5. NP = nasopharyngeal.

6. Does not distinguish between influenza A and B virus infections when used alone.

7. Immunofluorescence assay.

8. Requires use of a separate analyzer reader device

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1저자(연도)	기술문
Merckx (2017)	automated immunochromatographic antigen detection tests (digital immunoassays [DIAs]. Digital immunoassays use an instrument-based digital scan of the test strip to enhance antigen detection accuracy by eliminating the need for an operator to visualize and subjectively interpret test results.
Koski (2017)	the Sofia FIA uses an immunofluorescence-based lateral-flow technology to detect influenza A and B nucleoprotein antigens,
Chong (2019)	The Sofia Influenza A+B Fluorescence Immunoassay uses an immunofluorescence-based lateral flow technology to detect influenza A and B nucleoprotein antigen. This assay can be completed within 15 min and the results are reported as positive or negative for influenza virus A or B. Briefly, 250 µL of specimen and reagent buffer were mixed in a reaction tube. A final volume of 120 µL of sample was placed in a Sofia reaction cassette using a premeasured pipette and incubated for 15 min to allow the reaction to occur. The reaction cassette was inserted into the Sofia vfluorescent reader and the results were automatically printed within 1 min.
Yang (2018)	The Sofia influenza A+B fluorescent immunoassay (Sofia FIA) test utilizes immunofluorescence technology and a specific analyzing device to detect influenza nucleoproteins.
Ryu (2018)	The Sofia Influenza A+B Fluorescence Immunoassay (Quidel) is a rapid diagnostic kit that uses immunofluorescence technology to enhance its sensitivity . We compared three RIDTs with digital readout systems : BUDDI(NanoEnTek, Seoul, Korea), the Sofia Influenza A+B Fluorescence Immunoassay (Quidel, San Diego, CA, USA), and the BD Veritor System Flu A+B assay (BD Diagnostics, Sparks, MD, USA). A conventional rapid test (SD Bioline Influenza Ag A/B/A(H1N1/2009), Standard Diagnostics, Yongin, Korea) system for influenza detection was included for comparison.
Rack-Hoch (2017)	The new Sofia point-of-care system is a fluorescent immunoassay (FIA) based on lateral flow technology using monoclonal antibodies labelled with Europium as a fluorescent tag.
Ryu (2016)	Recently, influenza RDTs equipped with digital reader systems have

1저자(연도)	기술문
	been introduced and evaluated for comparison of detection sensitivity.[2,3] The Sofia® Influenza A + B Fluorescence Immunoassay (Sofia; Quidel Corp., San Diego, CA, USA) detects influenza A and B viruses using immunofluorescence technology with europium dye, and is interpreted with a portable fluorescence automatic reader.
Gomez (2016)	The Sofia Influenza A + B Fluorescence Immunoassay and SofiaRSV Fluorescence Immunoassay (FIAs) have recently been commercialised. They employ a novel technology different from other RDTs: these tests are immunofluorescence assays that use europium dye, which enhances the overall sensitivity. In addition, the assays have an instrument-read work mode, thus they avoid the operator interpretation of results that can lead to errors with other RDTs. The Sofia Influenza A + B and Sofia RSV FIAs (Quidel Corporation, San Diego, CA, USA) are lateral flow immunoassays , which employ a fluorescence reaction (europium dye) to detect influenzaA, influenza B or RSV nucleoprotein antigens.
Arbefeville (2015)	The Sofia Influenza A+B FIA is a novel fluorescence-based lateral flow immunoassay that uses specific antibody-coated beads coupled to a heat- and light-stable europium conjugate to detect Flu A and influenza B (Flu B) viral nucleoprotein antigens; this contrasts with the other assays, which are based on chromogenic or lateral-flow immunoassays. The Sofia Analyzer uses the europium fluorescent tag illuminated by an ultraviolet light source, to detect antigen. The fluorescence is measured by the Sofia Influenza Analyzer, which displays the test results as positive, negative, or invalid.
Noh (2015)	The Sofia™ Influenza A + B FIA (Quidel Corporation, CA, USA) is a novel fluorescent immunoassay used to detect influenza A and B within 15 minutes using the Sofia Analyzer
Leonardi (2015)	The Sofia influenza A B fluorescent immunoassay (FIA) (Quidel, Inc., San Diego, CA) uses a fluorescence reader to detect influenza nucleoprotein antigens.
Tuttle (2015)	First-generation rapid tests (laminar flow strip tests) produce macroscopic color changes visualizing the presence of viral antigens. The sensitivity of test strips is highly reader-dependent, and sensitivities vary widely between different study sites and products. In an effort to improve rapid test specificities and sensitivities, second-generation tests have been developed, utilizing advanced lateral flow technology, fluorescent tags, and automated readers for improved sensitivity and objectivity. => SOFIA
Hazelton (2015)	The Sofia Influenza A+B FIA (hereafter Sofia; Quidel, San Diego, CA) is a novel lateral-flow immunoassay that uses specific antibody-coated beads coupled to a heat and light-stable europium chelate instead of the traditional gold immunochromatography employed by most available RIDTs.
Bruning (2014)	Sofia Fluorescent Immunoassay Analyzer (FIA) (Quidel, San Diego, CA) is a rapid fluorescence-based lateral flow immunoassay in which results are analyzed by a compact instrument.
Dunn (2014)	The QuidelSofia® InfluenzaA+B FIA (Sofia; Quidel Corp., San Diego, CA, USA) detects and differentiates influenza A and B viruses using immunofluorescent technology to qualitatively identify virus-specific nucleoprotein antigens in nasal swabs, NP swabs, and NP aspirates or washes. Results from both the Veritor and Sofia assays are analyzed and interpreted digitallywith their respective readers. These digital immunoassays (DIAs) enable lower limits of detection (LODs) due to a low cut-off, reduce operator variability, and allow for a negative control that reduces the probability of non-specific binding. The objective, qualitative result of the DIAs potentially enables use of these assays outside the central laboratory environment.

1저자(연도)	기술문
Lewandrowski (2013)	Recently a novel fluorescent immunoassay (FIA) (Sofia Influenza A+B FIA, Quidel, San Diego, CA) was developed to detect and differentiate influenza A and B in patient specimens. The Sofia Influenza A+B FIA is a commercially available in vitro diagnostic device that detects viral nucleoprotein antigens. Although capable of detecting a broad range of influenza A virus subtypes, it was not designed to specifically identify the virus subtype. The test is a fluorescence-based lateral flow immunoassay in which results are analyzed objectively by a compact instrument (Sofia Analyzer) rather than visually. Instead of antibody-tagged colored microparticles, the Sofia technology uses a unique polystyrene microbead that has been impregnated with a chelate of europium.

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저자 (연도)	구분	A형		B형		A+B형		
		통합 Sn (95% CI)	통합 Sp (95% CI)	통합 Sn (95% CI)	통합 Sp (95% CI)	통합 Sn (95% CI)	통합 Sp (95% CI)	
Bruning (2017)	신속검사 전체	68.1 (58.0-76.0)	99.2 (95.5-99.6)	71.0 (56.8-82.1)	99.6 (99.2-99.8)	61.1 (53.3-68.3)	98.9 (98.4-99.3)	
	제품 별	Sofia	-	-	-	-	75.3 (59.2-91.5)	95.3 (91.5-99.2)
		QuickVue	-	-	-	-	44.6 (29.1-60.0)	99.3 (98.8-99.9)
		BinaxNow	-	-	-	-	44.1 (23.3-64.9)	99.4 (98.6-100)
		Directigen	-	-	-	-	35.8 (11.8-59.7)	99.2 (98.0-99.4)
		mariPOC	-	-	-	-	76.1 (53.5-98.7)	99.4 (98.3-100)
Merckx (2017)	traditional RIDTs	54.4 (48.9-59.8)	99.4 (99.1-99.7)	53.2 (41.7-64.4)	99.8 (99.7-99.9)	-	-	
	DIA	전체	80.0 (73.4-85.6)	98.3 (97.4-98.9)	76.8 (65.4-85.4)	98.7 (97.5-99.4)	-	-
		Sofia	77.8 (68.8-85.4)	98.5 (97.4-99.2)	73.5 (55.8-86.1)	98.0 (95.4-99.1)	-	-
		Veritor	83.0 (73.4-90.1)	97.5 (95.5-98.7)	80.0 (68.8-88.2)	99.5 (98.8-99.8)	-	-
	rapid NAATs	91.6 (84.9-95.9)	99.2 (98.6-99.7)	95.4 (87.3-98.7)	99.4 (98.9-99.8)	-	-	
Koski (2017)	Sofia (7편)	71.4-95.8	91.1-100	33.3-98.1	-	-	-	
	BD Veritor (3편)	72-93.8	57.1-99.07	69.3-94.2	100	-	-	
	Alere i* (6편)	73.2-99.3	62.5-100	45.2-100	53.6-100	-	-	

DIA, Digital immunoassay; RIDTs, Rapid Influenza Diagnostic Tests;

* iNAAT (isothermal nucleic acid amplification technology) 기술

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¶ Pandemic Influenza A(H1N1)pdm09 virus 검출

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- Green DA, George K. Rapid Antigen Tests for Influenza: Rationale and Significance of the FDA Reclassification. *Journal of Clinical Microbiology*. 2018;56: e00711-18.
- Hazelton BG, T.Ho, J.Ratnamohan, V. M.Dwyer, D. E.Kok, J. Detection of influenza A and B with the Alere TM i Influenza A & B: a novel isothermal nucleic acid amplification assay. *Influenza & Other Respiratory Viruses*. 2015a;9(3):151-4.
- Hazelton BN, G.Ratnamohan, V. M.Dwyer, D. E.Kok, J. Evaluation of the Sofia Influenza A + B fluorescent immunoassay for the rapid diagnosis of influenza A and B. *Journal of Medical Virology*. 2015b;87(1):35-8.
- Kammerer PER, J. M.Hawksworth, A. W.Myers, C. A.Brice, G. T. Performance of the Quidel Sofia rapid influenza diagnostic test during the 2012-2013 and 2013-2014 influenza seasons. *Influenza & Other Respiratory Viruses*. 2016;10(3):220-3.
- Koski RRR, M. E. A systematic review of rapid diagnostic tests for influenza: considerations for the community pharmacist. *Journal of the American Pharmacists Association: JAPhA*. 2017;57(1):13-9.
- Lee CKC, C. H.Woo, M. K.Nyeck, A. E.Lim, C. S.Kim, W. J. Evaluation of Sofia fluorescent immunoassay analyzer for influenza A/B virus. *Journal of Clinical Virology*. 2012;55(3):239-43.
- Leonardi GPW, A. M.Mitrache, I.Zuretti, A. R. Comparison of the Sofia and Veritor direct antigen detection assay systems for identification of influenza viruses from patient nasopharyngeal specimens. *Journal of Clinical Microbiology*. 2015;53(4):1345-7.
- Leonardi GPW, A. M.Zuretti, A. R. Comparison of conventional lateral-flow assays and a new fluorescent immunoassay to detect influenza viruses. *Journal of Virological Methods*. 2013;189(2):379-82.
- Lewandrowski KT, J.Menegus, M.Olivo, P. D.Lollar, R.Lee-Lewandrowski, E. Detection of influenza A and B viruses with the Sofia analyzer: a novel, rapid immunofluorescence-based in vitro diagnostic device. *American Journal of Clinical Pathology*. 2013;139(5):684-9.
- Merckx J, Wali R, Schiller I, Caya C, Gore GC, Chartrand C, Dendukuri N, Papenburg J. Diagnostic Accuracy of Novel and Traditional Rapid Tests for Influenza Infection Compared With Reverse Transcriptase Polymerase Chain Reaction. A Systematic Review and Meta-analysis. *Annals of Internal Medicine* 2017; 19:167(6):394-409.
- Noh JYC, W. S.Lee, J.Kim, H. L.Song, J. Y.Cheong, H. J.Kim, W. J. Clinical performance of the SofiaTM Influenza A+B FIA in adult patients with influenza-like illnesses. *Diagnostic Microbiology & Infectious Disease*. 2015;83(2):130-2.
- Olsen SJ, Kittikraisak W, Fernandez S, Suntarattiwong P, Chotpitayasunondh T. Challenges With New Rapid Influenza Diagnostic Tests. *Pediatr Infect Dis J*. 2014 January ; 33(1): 117-118.
- Piche-Renaud PPT, J.Chartrand, C.Gravel, J.Labrecque, M.Vallieres, E.Renaud, C. Evaluation of a fluorescent immunoassay rapid test (SofiaTM) for detection of influenza A+B and RSV in a tertiary pediatric setting. *Diagnostic Microbiology & Infectious Disease*. 2016;84(4):304-8.

- Rack-Hoch ALL, G.Hubner, J. Comparison of influenza and RSV diagnostic from nasopharyngeal swabs by rapid fluorescent immunoassay (Sofia system) and rapid bedside testing (BinaxNOW) vs. conventional fluorescent immunoassay in a German university children's hospital. *Infection*. 2017;45(4):529-32.
- Rath BT, F.Obermeier, P.Tuerk, E.Karsch, K.Muehlhans, S.Adamou, E.Duwe, S.Schweiger, B. Early detection of influenza A and B infection in infants and children using conventional and fluorescence-based rapid testing. *Journal of Clinical Virology*. 2012;55(4):329-33.
- Ryu SWS, I. B.Ryu, S. M.Shin, K. S.Kim, H. S.Kim, J.Uh, Y.Yoon, K. J.Lee, J. H. Comparison of three rapid influenza diagnostic tests with digital readout systems and one conventional rapid influenza diagnostic test. *Journal of Clinical Laboratory Analysis*. 2018;32(2).
- Ryu SWL, J. H.Kim, J.Jang, M. A.Nam, J. H.Byoun, M. S.Lim, C. S. Comparison of two new generation influenza rapid diagnostic tests with instrument-based digital readout systems for influenza virus detection. *British Journal of Biomedical Science*. 2016;73(3):115-20.
- Sakurai A, Shibasaki F. Updated values for molecular diagnosis for highly pathogenic avian influenza virus. *Viruses* 2012; 4: 1235-57.
- Selove WR, L. V. Performance of rapid SOFIA Influenza A+B test compared to Luminex x-TAG respiratory viral panel assay in the diagnosis of influenza A, B, and subtype H3. *Journal of Investigative Medicine*. 2016;64(4):905-7.
- Tuttle RW, A.Schwarz, W. S.Chen, X.Obermeier, P.Seeber, L.Tief, F.Muehlhans, S.Karsch, K.Peiser, C.Duwe, S.Schweiger, B.Rath, B. Evaluation of novel second-generation RSV and influenza rapid tests at the point of care. *Diagnostic Microbiology & Infectious Disease*. 2015;81(3):171-6.
- Yang JHH, P. Y.Shie, S. S.Yang, S.Tsao, K. C.Wu, T. L.Wu, T. S.Huang, C. T. Diagnostic performance of the Sofia influenza A+B fluorescent immunoassay in adult outpatients in Northern Taiwan. *Journal of Medical Virology*. 2018;90(6):1010-8.
- Yoon JY, S. G.Nam, J.Choi, S. H.Lim, C. S. The use of saliva specimens for detection of influenza A and B viruses by rapid influenza diagnostic tests. *Journal of Virological Methods*. 2017;243:15-9.

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검사원리 분류 신설 (2014-61호)	검사원리 분류 개정(2016.5.)	구체적인 검사법 예시
TIA	정밀면역검사	· TIA (Turbidimetric immunoassay)
		· NIA (Nephelometric immunoassay) 등
		· EIA (Enzyme immunoassay)
		· ELISA (Enzyme linked immunosorbent assay)
		· EMIT (Enzyme multiplied immunoassay technique)
		· ELISPOT (enzyme-linked immunosorbent spot)
EIA		· ELIFA (Enzyme linked immunofiltration assay)
		· ELFA (Enzyme linked fluorescent assay)
		· MEIA (Microparticle enzyme immunoassay)
		· Line immunoassay (Immunoblot)
		· FEIA (Fluorescence enzyme immunoassay)
		· SBPA (Sandwich binding protein assay)
		· LIFA (Ligand-mediated immunofunctional assay)
		· Immuno-filtration assay 등
CIA		· CIA (Chemiluminescence immunoassay)
		· CMIA (Chemiluminescence microparticle immunoassay)
		· ECLIA (Electrochemiluminescence immunoassay) 등
FIA		· FIA (Fluorescence immunoassay)
		· FMIA (Fluorescence microbead immunoassay)
		· FPIA (Fluorescence polarization immunoassay) 등
RIA		· RIA (Radioimmunoassay)
		· IRMA (Immunoradiometric assay)
		· RRA (Radio receptor assay) 등
Neutralization test	Neutralization test	Neutralization test
EP	EP	· EP (Electrophoresis)
Column chromatography	정밀분광/질량분석	· GC (Gas chromatography)
		· LC (Liquid chromatography)
		· HPLC (High performance liquid chromatography) 등
MS		· MS (Mass spectrometry)
AAS		· AAS (Atomic absorption spectroscopy)
		· AAM (Atomic absorption method) 등
ICP		· ICP-ES/-MS (Inductively coupled plasma-Emission spectroscopy/-mass spectrometry) 등
Extended chromatography		· Chromatography-MS
		· MS-MS
		· Chromatography-MS-MS 등
Western blot	Western blot	· Western blot 등

검사원리 분류 신설 (2014-61호)	검사원리 분류 개정(2016.5.)	구체적인 검사법 예시
Immunofluorescent assay (IFA)	Immunofluoresce nt assay (IFA)	<ul style="list-style-type: none"> · DIFA (Direct immunofluorescent assay) · IIFA (Indirect immunofluorescent assay) · Fluorescent antibody test · FTA-ABS · FANA (Fluorescent antinuclear antibody assay) 등
IEP	면역전기영동	· IEP (Immunoelectrophoresis)
IFE		· Rocket EP (Laurell method) 등
Flow cytometry	Flow cytometry	· Flow cytometry 등

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